

BUGBOX SPECIFICATION SUMMARY

| MODEL | | BUGBOX/BUGBOX M | BUGBOX PLUS |
|-------------------------------|--------------|------------------|-------------|
| External Dimensions | Width | 800 mm | 830 mm |
| | Depth | 660 mm | 660 mm |
| | Height | 650 mm | 650 mm |
| Internal Dimensions | Width | 500 mm | 500 mm |
| | Depth | 460 mm | 460 mm |
| | Height | 420 mm | 420 mm |
| Maximum Capacity | 90 mm Plates | 270 | 234 |
| Working Capacity | 90 mm Plates | 200 | 180 |
| Interlock Dimensions | Width | 100mm | 150 mm |
| | Depth | 100mm | 230 mm |
| | Height | 200mm | 190 mm |
| Interlock Capacity | 90 mm Plates | 10 | 18 |
| Interlock Time Cycle | | 15 sec | 35 sec |
| Interlock Door Operation | | Manual | Manual |
| Weight* | | 99 lbs / 143 lbs | 121 lbs |
| Petri Dish Holders (Standard) | | 3 | 3 |


*Bugbox M (143lbs) includes ICONIC™, the gas mixing system from Baker Ruskinn that gives Bugbox M unparalleled oxygen control while using up to 40% less nitrogen.

STANDARD FEATURES

- Detox advanced carbon filtration system
- Ezee Sleeve™ direct hand entry system
- Energy saving fluorescent illumination
- Inspection spot lamp
- Low gas alarm
- Automatic humidity control
- Palladium catalyst
- Anaerobic indicator strips
- Petri dish holders

OPTIONS & ACCESSORIES

- Vacuum line port
- Gas sample port
- Cable gland port
- Internal electrical outlet
- Gas tank regulators and filter modules
- Workstation stand
- External docking facility for anaerobic jars
- Power failure back-up system
- Data logging connection
- Single Plate Entry System (SPES)
- For facultative and microaerophiles Bugbox-M allows user defined control of O₂ and CO₂



ICONIC™ allows:

- O₂ from 0.0% to 23.0% in 0.1% increments
- CO₂ control from 0.1% to 30.0% in 0.1% increments
- O₂ sensor calibration with one touch
- Microaerophilic cycling, facilitating up to 4 different O₂ and CO₂ concentrations through a user-defined sequence of time

Concept

Anaerobic & Microaerophilic Workstations

If you are looking for more robust processing power and capacity, the Baker Ruskinn Concept range of workstations are the perfect addition to any lab.



Scan the QR code above to read more

BAKER RUSKINN

A NEW LEVEL OF VERSATILITY, EFFICIENCY AND POWER IN A COMPACT SPACE



Get in touch today to hear how we can work with you.

For U.K. / global inquiries (except U.S./Canada):
sales@ruskinn.com
+44 (0) 1656 645988
www.bakerco.com

Connect and Collaborate with Baker



Designed in the U.K.

Bugbox

ANAEROBIC & MICROAEROPHILIC WORKSTATIONS
GLOBAL BROCHURE (220;110V)

SEE THINGS DIFFERENTLY

Experience Baker Ruskinn Anaerobic Workstations



VERSATILE AND FLEXIBLE TO FIT YOUR WORKLOAD

Multiple models and a variety of options are available to fit your specific needs.

BUGBOX - YOUR PERSONAL WORKSTATION

- Up to 270 90mm plate capacity
- Interlock transfer - 10 plates in 15 seconds

BUGBOX PLUS - COMPACT WORKSTATION WITH LARGER INTERLOCK

- Up to 234 90mm plate capacity
- Interlock transfer - 18 plates in 35 seconds

BUGBOX M

- For facultative and microaerophiles
- User defined control of O₂ featuring ICONIC™ from Baker Ruskin

Bugbox anaerobic workstations are designed to help microbiologists cope with rising workloads and provide the best primary isolation rates. Plates can be examined easily without exposing them to oxygen. The interlock system allows simple and fast transfer of 90mm plates into the anaerobic chamber.

With quick and easy access via the Ezee Sleeve™ Glove port system and energy-saving lighting that provides perfect illumination, Bugbox is easy to use. Its compact size meets the needs of even the smallest laboratory spaces. Adjustable temperature and humidity provides a precisely controlled anaerobic environment that is optimal for cell growth, with no dry spots.

Compared with approximately 20 anaerobic jars per week, Bugbox is economical with a lower cost per plate, more reliable, providing a stable atmosphere, and minimal maintenance.

For 20 years Baker Ruskinn anaerobic workstations have been the trusted choice for laboratories around the world. More than 1000 anaerobic workstations are installed in more than 40 countries - and more than 200 research publications feature the Baker Ruskinn anaerobic technology.

Our anaerobic workstations are designed to help microbiologists cope with rising workloads and provide the best primary isolation rates.

DESIGNED TO PROTECT YOUR RESULTS

The acrylic airtight chamber is flooded with anaerobic gas mix (H₂ in N₂) and O₂ is displaced.

If any O₂ remains or is allowed to enter, it is "scavenged" by a palladium catalyst situated under the floor tray - the O₂ reacts with the H₂ to form water.

Interlock uses an N₂ purge, so when a user brings in plates through the interlock, no O₂ enters the main chamber - inner and outer interlock doors cannot be opened simultaneously.

Gloveless Ezee Sleeves™ are purged using N₂ gas via foot pedals, so no O₂ enters the main chamber when the glove ports are opened.

CONVENIENT & COMFORTABLE USER EXPERIENCE

1 QUICK AND EASY DIRECT ACCESS

Gloveless, cuffed sleeve system (Ezee Sleeve™) takes less than 40 seconds for direct hand access into the chamber.

2 SHORTEST INTERLOCK CYCLE TIME IN THE INDUSTRY

As little as 15 seconds for the 10 plate capacity interlock.

3 SINGLE PLATE ENTRY SYSTEM (SPES)

This optional accessory is a mailbox like slot, which allows quick side entry or exit of individual plates, bypassing the interlock cycling process.

ECONOMIC AND RELIABLE FOR LONG TERM SAVINGS

- Standard dual gas operation, low gas consumption and running costs.
- Lower cost per plate compared to anaerobic jars.
- Minimal maintenance and downtime.

ULTIMATE CONTROL FOR OPTIMUM CELL ENVIRONMENT

- Accurate temperature control from ambient + 5°C to 45°C.
- Accurate and automated humidity control, no dry spots.
- Palladium catalyst maintains anaerobic environment, plus anaerobic colour-indicator strips verify anoxic conditions.
- Ezee Sleeve™ Direct Hand entry system allows access without disrupting the atmosphere within the chamber.

