

EXTERNAL CONDITIONING UNITS

Sometimes due to geometry, size or the elements that compose the system to be tested, is difficult to introduce the equipment within a standard climatic chamber.

For testing of climate conditioning in inaccessible parts are often used external generators, which have associated an external workspace of small dimensions, that allow conditioning this particular area of the piece.

Generator system can be built to acclimate the box in temperature or temperature and humidity. And tests box is manufactured adapted to the system to be tested.

Compact chambers can be modified so that they function as generator system, using universal port placed in the side and an adaptable system to connect recirculation hoses. The chamber could work like a chamber or like a generator.



- Conditioning in temperature or in temperature and humidity.
- Custom made external box.
- Different thermal gradient options.
- Insulating hoses with different lengths.
- In accordance with customer's need.