



KOJAIR®

# CLEANWIZARD H

## Product Safety Horizontal flow Cabinet

### GENERAL

Safe, reliable and economic safety cabinet for protection of the work item including ergonomic practical properties needed in normal work.

### TECHNICAL INFORMATION

#### Materials and finishing

All materials and surfaces are smooth and resistant to cleaning and disinfection agents.

#### Lighting

Glare-free and shadow-free fluorescent lighting providing illumination on work surface 1200 lux

#### Standard equipment

- Flow control with digital display
- Alarms for false flow speed
- 1/2 speed via switch
- UV outlet via switch
- Base stand 700, 800 or 900 mm



#### Filters

HEPA H 14 99,995 % MPPS  
acc. to EN 1822

#### Air speed

Air flow 0,4 m/s ± 20%

#### Sound level

The sound level with clean filters is 55 dB.

#### Options

- UV light 30 W
- Gas tap
- Socket

CleanWizard H	130 H	200 H
Total width ( mm )	1280	1910
Depth	870	870
Working space ( mm )		
Width	1180	1830
Depth	500	500
Height	610	610
Power ( max. W )	470	730
Weight ( incl. base stand kg )	180	220

**Manufacturing Plant and Head Office:** Teollisuustie 3, FIN-35700 Vilppula, Finland

Tel. + 358-(0)3-471 7000,

Fax + 358-(0)3-471 7077, <http://www.kojair.com>, [sales@kojair.com](mailto:sales@kojair.com)