



Save your brain, concentrate on work

BIOWIZARD Platinum NRS

Microbiological SMART Safety Cabinet





EN 12469:2000 Norm

Current Microbiological Safety Cabinets have high noise levels and the EN norm settles for under 65dB. Noise is undesired and unpleasant for workers. It puts load on brains and can harm physically and mentally. High sound level not only affects hearing, stress level and blood pressure but also weakens concentration, creativity and comprehension. The Kojair Platinum SMART is build conform the EN 12469:2000, with a sound level that starts at 37dB!



Noise reduction system (NRS)

Thanks to the Kojair noise reduction system (NRS) the Kojair Platinum Smart microbiological safety cabinet has a sound level from 37dB. This is based on the functional and safety principles of a microbiological safety cabinet. Safety comes first!!



Officially tested and certified

More and more manufacturers claim to have microbiological safety cabinet with a low noise level. But after purchase the microbiological safety cabinet doesn't reach the indicated values. To prevent this discussion Kojair's Platinum SMART cabinet has been officially tested and certified by an independent agent, Technical Research Centre of Finland Ltd (VTT). Tests have been executed conform to the EN 12469 / NEN-EN-ISO 3744.







