

MF SERIES

VERTICAL PARALLEL LIFT DOOR · REFINED DESIGN · CERAMIC FIBER INSULATION (1200/1400, 300/400)

High Quality 24/7 Continuous-T° Muffle Furnaces · from 1000°C to 1400°C

STANDARD FEATURES

- CE manufactured
- Maximum operating temperature: 1000°C - 1200°C - 1300°C - 1400°C
- Rapid heating
- Compact and lightweight
- 24/7 continuous work capability
- 4 heating plates with KANTHAL AF and APM
- Low thermal mass insulation
- Built with low density ceramic bricks and ceramic fiber
- Double insulation includes air chamber
- Outer case in painted metal sheet (inox optional)
- Rear ventilation via ceramic chimney
- Ceramic tray included
- Thermocouple type K or S
- Spare parts easily replaceable by end user

FURNACE CONTROLS

- Lower front control panel
- Solid state relay
- General safety switch
- General safety contactor

PAD Digital control

- PID parameters
- Non-volatile memory
- Microprocessor-based temperature controls
- Alarm

CONTROL OPTIONS

- Ramp programmers up to 64 segments
- Optional:
- Eurotherm EPC Series + Data logger Itools sofware by Ethernet
- Eurotherm Nanodac Series +Data logger and Itools sofware by Ethernet according AMS2750E and 21CFR Part 11

SAFETY SHUT-OFF

- Thermocouple break shut-off
- Turns off upon door opening

ACCESSORIES

- Interchangeable temperature-uniform trays with rim
- Refractory, incoloy, SS and ceramic tray
- Forced cooling system
- Smoke chimney
- Forced smoke chimney
- Safety alarm Class II. Over-temperature protection
- Inlet gas entry
- Flow meter box
- and more, ask for our full assortment!





MF SERIES

VERTICAL PARALLEL LIFT DOOR · REFINED DESIGN · CERAMIC FIBER INSULATION (1200/1400, 300/400)

High Quality 24/7 Continuous-T° Muffle Furnaces · From 1000°C to 1400°C

CHARACTERISTICS

- Modern design metal case with chrome-phosphatizing base protection and external finish with heat-resistant metal paint.
- Refractory parts engineered to resist extreme temperature changes, and specific ceramic paste types applied according to temperature and work fatigue of each part.
- Heat resistance in refractory insulation of very low thermal conductivity coefficient.
- Door system adjusted on the furnace frame by pressure, allowing for complete sealing. Electrically and thermally insulated door handle.

ACCESSORIES

- Extraction Chimney: Self-extraction design to eliminate smoke in processes that produce smoke in a considerable amount or when smoke extraction is advisable due to the nature of the process. Chimney outlet connection to a smoke bell or to the exterior by end user.
- Forced air extraction chimney: Specially designed for a forced self-extraction to evacuate smoke fastat resistance in refractory insulation of very low thermal conductivity coefficient.
- · Bottom trays: Interchangeable, temperature uniform, with rim to protect against spilling, fusion or adherence of materials.

SPECIFICATIONS

Fully customized solutions by request We reserve the right to change technical specifications

	Inner dimensions Outer dimensions						Maximum	Maximum	Maximum	Termo-	Control	Heating				
Model	mm		mm			Volume	Power	Voltage	Heated	Temperature	Temperature ° C	Temperature	couple	Type	Elements	
	High	Wide	Deep	High	Wide	Deep	Liters	Kw	v	Zones	۰c	limited time	Continuous °C			
MF3-124	150	150	200	620	430	470	3	2,5	220	4	1200	1150	1100	K	R. Prog	Kanthal Wire
MF6-124	150	150	250	530	470	490	6	2,5	220	4	1200	1150	1100	K	R. Prog	Kanthal Wire
MF9-124	150	200	300	680	470	600	9	5,5	220	4	1200	1150	1100	K	R. Prog	Kanthal Wire
MF12-124	150	200	400	680	470	730	12	5,5	220	4	1200	1150	1100	K	R. Prog	Kanthal Wire
MF22-124	200	280	400	730	550	750	22	8,8	380 III	4	1200	1150	1100	K	R. Prog	Kanthal Wire
MF30-124	250	280	400	770	550	750	30	8,8	380 III	4	1200	1150	1100	K	R. Prog	Kanthal Wire
MF40-124	250	280	500	770	550	850	40	12,5	380 III	4	1200	1150	1100	K	R. Prog	Kanthal Wire

Inn		Inner dimensions		Outer dimensions							Maximum	Maximum	Maximum	Termo-	Control	Heating
Model	mm		mm			Volume	Power	Voltage	Heated	Temperature	Temperature ° C	Temperature	couple	Type	Elements	
	High	Wide	Deep	High	Wide	Deep	Liters	Kw	v	Zones	۰c	limited time	Continuous °C			
MF3-134	150	150	200	620	430	470	3	2,5	220	4	1300	1250	1200	S	R. Prog	Kanthal Wire
MF6-134	150	150	250	530	470	490	6	2,5	220	4	1300	1250	1200	S	R. Prog	Kanthal Wire
MF9-134	150	200	300	680	470	600	9	5,5	220	4	1300	1250	1200	S	R. Prog	Kanthal Wire
MF12-134	150	200	400	680	470	730	12	5,5	220	4	1300	1250	1200	S	R. Prog	Kanthal Wire
MF22-134	200	280	400	730	550	750	22	8,8	380 III	4	1300	1250	1200	S	R. Prog	Kanthal Wire
MF30-134	250	280	400	770	550	750	30	8,8	380 III	4	1300	1250	1200	S	R. Prog	Kanthal Wire
MF40-134	250	280	500	770	550	850	40	12,5	380 III	4	1300	1250	1200	S	R. Prog	Kanthal Wire

	Inner dimensions			Outer dimensions							Maximum	Maximum	Maximum	Termo-	Control	Heating
Model	mm		mm			Volume	Power	Voltage	Heated	Temperature	Temperature ° C	Temperature	couple	Type	Elements	
	High	Wide	Deep	High	Wide	Deep	Liters	Kw	v	Zones	۰c	limited time	Continuous °C			
MF3-144	150	150	200	620	430	470	3	2,5	220	4	1400	1350	1300	S	R. Prog	Kanthal Wire
MF6-144	150	150	250	530	470	490	6	2,5	220	4	1400	1350	1300	S	R. Prog	Kanthal Wire
MF9-144	150	200	300	680	470	600	9	5,5	220	4	1400	1350	1300	S	R. Prog	Kanthal Wire
MF12-144	150	200	400	680	470	730	12	5,5	220	4	1400	1350	1300	S	R. Prog	Kanthal Wire
MF22-144	200	280	400	730	550	750	22	8,8	380 III	4	1400	1350	1300	S	R. Prog	Kanthal Wire
MF30-144	250	280	400	770	550	750	30	8,8	380 III	4	1400	1350	1300	S	R. Prog	Kanthal Wire
MF40-144	250	280	500	770	550	850	40	12,5	380 III	4	1400	1350	1300	S	R. Prog	Kanthal Wire

Up to 1000°C available

