



SCI-tive Dual Specifications

R = Right-hand chamber

L= Left-hand chamber

	SCI-tive Dual Symmetrical	SCI-tive Dual Asymmetrical
Workstation Fabrication	Solvent Bonded Acrylic	
External Dimensions	Width: 2995 mm / 111"	Width: 2815 mm / 118"
	Depth: 826 mm / 33"	Depth: 975 mm / 38"
	Height (excluding stand): 1077 mm / 42"	Height (excluding stand): 1077mm – 1195 mm / 42" – 47"
Internal Dimensions	Width: 1200 mm / 47" (L & R)	Width: 1200 mm / 47" (L) and 1020 mm / 40" (R)
	Depth: 555 mm / 24" (L & R)	Depth: 555 mm / 24" (L) and 800 mm / 30" (R)
	Height: 690 mm / 27" (L & R)	Height: 690 mm / 27" (L) and 808 mm / 32" (R)
Workstation Weight	approx: 355 kg / 783 lbs.	
Interlock Dimensions	Width: 270 mm / 11"	
	Depth: 200 mm / 8"	
	Height: 210 mm / 8"	
Interlock Cycle Time	60 seconds	