


Skintact® Electrosurgical Neutral Electrode

RSW23A30

Code: 50907



PRODUCT	
Characteristics	disposable, self adhesive, non-sterile
Type	split
Shelf Life	36 Mon
Storing Conditions (min./max.)	+5°C (41°F) / +30°C (86°F)
Patient Target group	Neonatal (0.45 - 4.99 kg / 1 - 11 lbs)
Cable length [m]	3
Connector	REM

CLASSIFICATION AND STANDARDS	
Classification (REG(EU) 2017/745, Annex VIII)	Class IIb

BIOCOMPATIBILITY	
ISO 10993-1 ff	pass

PERFORMANCE	
IEC 60601-2-2	pass

MATERIALS	
Carrier Material	Non woven Medical Tape
Adhesive	Polyacrylate
Conductive Layer	Aluminium
Gel	Hydrogel
Protective Cover	Siliconized Foil

DIMENSIONS	
Electrode Size L/W [cm] (approx.)	17 / 8
Adhesive Area [cm²] (approx.)	86
Conductive Area [cm²] (approx.)	43
Gel Thickness [mm]	0.5 ± 0.05

STANDARD PACKAGING	
Pieces / Pouch	1
Pieces / Box	10
Pieces / Case	40
Duplex Label on Pouch	yes
Lot number and expiration date on every pouch, box and case.	

DIMENSIONS PACKAGING	
Pouch (L/W) [mm]	260 x 196
Box (L/W/H) [mm]	282 / 185 / 141
Case (L/W/H) [mm]	398 / 298 / 304

MATERIAL PACKAGING	
Pouch	Paper / Aluminium / PE
Box	Cardboard
Case (L/W/H) [mm]	Cardboard

ENVIRONMENTAL INFORMATION	
Compliance with RoHS Directive	yes
Compliance with REACH Regulation	yes
Not made with natural rubber latex.	

POWER & CURRENT INFORMATION	
Determine Watt maximum for each mode of operation as stipulated in IFU of ESU device manufacturer (maximum current at 350 mA)	