


Skintact® Electrosurgical Neutral Electrode

WR21A30

Code: 59597



PRODUCT

Characteristics	disposabel, self adhesive, non-sterile
Type	split
Shelf Life	36 Mon
Storing Conditions (min./max.)	+5°C (41°F) / +30°C (86°F)
Patient Target group	Universal (> 5 kg / > 11 lbs)
Cable length [m]	3
Connector	REM

STANDARD PACKAGING

Pieces / Pouch	1
Pieces / Box	50
Pieces / Case	200
Duplex Label on Pouch	yes
Lot number and expiration date on every pouch, box and case.	

CLASSIFICATION AND STANDARDS

Classification (REG(EU) 2017/745, Annex VIII)	Class IIb
Classification (21 CFR 878.4400)	Class II
CMDR SOR/98-282	Class II

DIMENSIONS PACKAGING

Pouch (L/W) [mm]	210 / 196
Box (L/W/H) [mm]	385 / 285 / 132
Case (L/W/H) [mm]	400 / 300 / 570

BIOCOMPATIBILITY

ISO 10993-1 ff	pass
----------------	------

PERFORMANCE

IEC 60601-2-2	pass
---------------	------

MATERIAL PACKAGING

Pouch	Paper / Aluminium / PE
Box	Cardboard
Case	Cardboard

MATERIALS

Carrier Material	Non woven Medical Tape
Adhesive	Polyacrylate
Conductive Layer	Aluminium
Gel	Hydrogel
Protective Cover	Siliconized Foil

ENVIRONMENTAL INFORMATION

Compliance with RoHS Directive	yes
Compliance with REACH Regulation	yes
Not made with natural rubber latex.	

DIMENSIONS

Electrode Size L/W [cm] (approx.)	16 / 13
Adhesive Area [cm²] (approx.)	163
Conductive Area [cm²] (approx.)	90
Gel Thickness [mm]	0.5 ± 0.05