


#### Skintact® Electrosurgical Neutral Electrode

**WR21A50**

Code: 59980



#### PRODUCT

Characteristics	disposabel, self adhesive, non-sterile
Type	split
Shelf Life	36 Mon
Storing Conditions (min./max.)	+5°C (41°F) / +30°C (86°F)
Patient Target group	Universal (> 5 kg / > 11 lbs)
Cable length [m]	5
Connector	REM

#### STANDARD PACKAGING

Pieces / Pouch	1
Pieces / Box	50
Pieces / Case	200
Duplex Label on Pouch	yes
Lot number and expiration date on every pouch, box and case.	

#### CLASSIFICATION AND STANDARDS

Classification (REG(EU) 2017/745, Annex VIII)	Class IIb
Classification (21 CFR 878.4400)	Class II
CMDR SOR/98-282	Class II

#### DIMENSIONS PACKAGING

Pouch (L/W) [mm]	210 / 196.5
Box (L/W/H) [mm]	385 / 285 / 160
Case (L/W/H) [mm]	395 / 595 / 358

#### BIOCOMPATIBILITY

ISO 10993-1 ff	pass
----------------	------

#### MATERIAL PACKAGING

Pouch	Paper / Aluminium / PE
Box	Cardboard
Case	Cardboard

#### PERFORMANCE

IEC 60601-2-2	pass
---------------	------

#### ENVIRONMENTAL INFORMATION

Compliance with RoHS Directive	yes
Compliance with REACH Regulation	yes
Not made with natural rubber latex.	

#### MATERIALS

Carrier Material	Non woven Medical Tape
Adhesive	Polyacrylate
Conductive Layer	Aluminium
Gel	Hydrogel
Protective Cover	Siliconized Foil

#### DIMENSIONS

Electrode Size L/W [cm] (approx.)	16 / 13
Adhesive Area [cm²] (approx.)	163
Conductive Area [cm²] (approx.)	90
Gel Thickness [mm]	0.5 ± 0.05