




QUALITY THAT WEIGHS





The entry level of balances keeping all the quality of higher level balances. This model has all necessary features to obtain precise and reliable results with a low price. The refined design, the small dimensions, the backlit LCD display, the possibility of operating also with AA type batteries (not included) and the many on-board software functions already included make this balance a product with the best quality over price.

Main features

- High readability load cell weighing system
- User selectable environmental filter
- Small footprint, portable
- External calibration
- Backlight LCD display
- Operates also with AA batteries (not included)
- Plastic protection cover

Accessories

- Factory Calibration certificate (BL0333)

Technical data

- Span drift (+ 10...+ 30 °C): +/- 6ppm/°C
- Dimensions LxWxH (mm): 190x255x80
- Power supply 110-230Vac, 50/60Hz; output 9V 1A multi plug
- Net Weight: 1 Kg



Model	Max (g)	Readability [d] (g)	Pan (mm)	Repeatability (St.Dev.) (g)	Linearity (g)	Response time (sec.)
-------	---------	---------------------	----------	-----------------------------	---------------	----------------------

External calibration

ES422	420	0,01	Ø 110	0,01	± 0,02	≤ 2
ES2201	2200	0,1	150x140	0,1	± 0,3	≤ 2