



thermostatic baths

WWW.FALCINSTRUMENTS.IT

water baths

HEATING IMMERSION UNITS



The temperature is set and checked by a digital thermostat with PID controller and PT100 probe. The stirring pump ensures excellent temperature uniformity inside the tank



The instrument is equipped with a safety float that locks the adjustment in the event that the liquid level is too low



Clamp included



FA 90
with internal pump

Technical data		FA 90
Pump capacity	l/min	9
Heating power	W	1300
Precision	°C	± 0,5
Head	m	1,8
Temperature	°C	+5 amb. to +150
Digital timer	hh:mm	∞ or 99:59
Display resolution	°C	0.1
Immersed pump height	mm	120
External dimensions LxWxH	mm	160x140x270
Weight	kg	3
Code		612.1029.00

V/HZ 220-240/50-60

SAFETY CLASS 3.1

PROTECTION CLASS IP 30

ACCESSORIES



EXTERNAL PUMP SET
It enables the external power supply via a closed circuit. Ø12mm tube.
Code 612.1029.02



18/8 stainless steel tank
with bridge

Dimensions	lt	Code
330x180x180	5	612.1029.26
360x350x180	12	612.1029.29
530x350x150	18.5	612.1029.32



Plexiglass tank
max oper. temp. 60°C

Dimensions	lt	Code
120x410x150	7	612.1029.35
300x330x150	12	612.1029.37
300x500x150	18.5	612.1029.39



18/8 stainless steel
lid with handle

Tank lt	Code
5	612.1029.27
12	612.1029.30
18.5	612.1029.33