



mini ovens & mini incubators

WWW.FALCINSTRUMENTS.IT

mini incubators

MINI INCUBATORS



ICT



Incubators with transparent Plexiglas door to inspect the chambers



Digital temperature control via PID action with PT 100 probe



Natural ventilation



Ø 10 mm hole on the top side to insert a thermometer or a probe

Technical data		ICT 5,4	ICT 18	ICT 28
Capacity	l	5,4	18,0	28,0
Min temperature	°C	room+ 5	room+ 5	room+ 5
Max Temperature	°C	+80	+80	+80
Precision	37°C	±0,5	±0,5	±0,5
Resolution	°C	0,1	0,1	0,1
Heating power	W	120	360	420
Included Shelves	No.	2	2	2
Max no. of shelves	No.	4	5	5
Internal dimensions LxWxH	mm	190x150x190	300x300x200	400x250x300
External dimensions LxWxH	mm	400x260x300	520x390x320	640x360x440
Weight	kg	10	15	26,5
Code		720.2070.10	720.2070.12	720.2070.14

V/HZ 220/240-50/60

SAFETY CLASS 0

PROTECTION CLASS IP 42