



mini ovens & mini incubators

WWW.FALCINSTRUMENTS.IT

mini incubators

MINI PORTABLE INCUBATORS WITH COOLING AND HEATING



Ideal for measurements and tests on the field
(also for bench applications in laboratory)



Digital timer: option to program with 5
temperature steps
Cooling via Peltier technology



Power supply 220V and 12V (cigarette lighter)



ICT-P

Technical data		ICT-P
Capacity	l	18
Temperature	°C	10 to 65
Precision	37°C	±0,2 °C
Resolution	°C	0,1
Uniformity	°C	0,1
Stability	°C	±0,2
Timer		- : - - (∞)
	mm:ss	00:01 to 59:59
	hh:mm	1:00 to 23:59
	dd:hh	1:00 to 45:00
Timer resolution	s	1
Heating and cooling setting		automatic detection
Internal dimensions LxWxH	mm	275x185x335
External dimensions LxWxH	mm	295x350x465
Weight	kg	6
Power Supply	W	60
Power	V	12
Code		720.2070.00

V/HZ 100 TO 220/240/50-60

SAFETY CLASS 0

PROTECTION CLASS 00



FALC INSTRUMENTS s.r.l.

Via G. M. Compagnoni, 2
24047 Treviglio (BG) - Italy
+39 0363 304660 | falc@falcinstruments.it

WWW.FALCINSTRUMENTS.IT