



# thermostatic baths

[WWW.FALCINSTRUMENTS.IT](http://WWW.FALCINSTRUMENTS.IT)

# refrigerated baths

## CIRCULATION-COOLING BATHS



SBF



LCD touch screen display  
Tanks are made of one piece  
of 18/8 stainless steel sheet



The temperature is checked by a PID microprocessor  
and PT100 probe



Automatic shutdown of the compressor over 37.0°C



Digital timer hh:mm:ss or infinite



Eco Friendly function, the cooling compressor starts working  
only when necessary



Easy access to keep the instrument clean and effective at all times  
Manual operation of the fireplace  
Setting for 4 languages (Italian-English-French-Spanish)  
Program saving function (temperature/maintenance/pump speed)  
Programmable delayed start (day hour)

### Technical data

|                         |          | SBF 7T             | SBF 15T            |
|-------------------------|----------|--------------------|--------------------|
| Capacity                | lt       | 7                  | 15                 |
| Drain tap               |          | no                 | no                 |
| Timer                   | hh:mm:ss | ∞ or 99:59:59      | ∞ or 99:59:59      |
| Working temperature     | °C       | -30 to +150        | -30 to +150        |
| Precision in the tank   | °C       | ±0.2               | ±0.2               |
| Heating capacity        | W        | 1300               | 1300               |
| Cooling capacity        | °C       | -25 0 25           | -25 0 25           |
|                         | W        | 210 170 260        | 520 600 850        |
| Pump type               |          | Pressure           | Pressure           |
| Pump Flow rate          | lt/min   | 7.8                | 4                  |
| Prevalence              | mt       | 2.5                | 2.5                |
| Bottom dimensions LxWxH | mm       | 100x250x200        | 250x290x200        |
| Usable tank size LxWxH  | mm       | 100x150x200        | 250x190x200        |
| Outer dimensions LxWxH  | mm       | 300x440x730        | 440x475x730        |
| Weight                  | kg       | 31                 | 41                 |
| <b>Code</b>             |          | <b>612.1050.08</b> | <b>612.1050.16</b> |

V/HZ 220-240/50-60

SAFETY CLASS 2

PROTECTION CLASS IP 30

### ACCESSORIES



Stainless steel floating level device for water  
shortage

**Code 612.1029.04**

Wheels

**Code 612.1050.31**

External circulation pipes 1.5 m

**Code 612.1050.33**

1 lt liquid for cryostat -30÷80°C 1 LT

**Code 612.1050.35**