



thermostatic baths

WWW.FALCINSTRUMENTS.IT

water baths

THERMOSTATIC BATHS WITH DUBNOFF TYPE SHAKING ACTION



WB-MF



Ideal for the control of samples, incubation, Cell Cultivation, Food and Beverage Testing



The temperature is set and checked by a digital thermostat with PID controller and PT100 probe



Metal structure with anticorrosive coating. The tanks are made of 18/8 stainless steel sheets pressed in a single piece



Digital speed control



Manual reset safety thermostat

Technical data

		WB-MF
Capacity	lt	20
Oscillation amplitude	mm	15
Oscillation frequency	Rpm	0-200
Heating capacity	W	2100
Precision	°C	±0.5
Temperature	°C	+5 amb. to +200
External dimensions LxWxH	mm	670x400x290
Internal dimensions LxWxH	mm	440x240x200
Weight	kg	16
Code		610.1007.74

V/HZ 220-240/50-60

SAFETY CLASS 1

PROTECTION CLASS IP 42

ACCESSORIES



Polycarbonate gable lid with hinge

Code 610.1007.76



General purpose container
WxD 325x235x60 mm
Code 610.1007.81



Supporting plate without clips
for flasks and beakers

Capacity	LxP 365x235 mm
ml	max holder pieces for plates
25 e 50	24-20
100	15
250	12
500	6
1000	4
2000	2
Code	610.1007.83



Fixing clips for flasks
and beakers

Capacity	Code
ml 100	603.0131.75
ml 250	603.0131.76
ml 500/1000	603.0131.77
ml 2000	603.0131.78