

Circulaire® Touchscreen Fume Cupboards



CT 1100

CT PRO 1400

**BRITISH
MANUFACTURED**



**FIVE YEAR
WARRANTY**



**ULTRA DEEP
CARBON
FILTERS**



**CUSTOMISABLE
FOR YOUR
APPLICATION**



The Circulaire® Touchscreen Recirculating Filtration Fume Cupboards feature the very latest in air filtration technology.

Fully compliant with COSHH regulations and international standards, the cupboards are certified to BS7989:2001 for Filtration Fume/Particulate Cupboards.

The units employ extra-large capacity Activated Carbon filters for removal of fumes and/or HEPA filters for particulate removal. The filters are up to 30% bigger than those of rival Fume Cupboards and are significantly deeper making fume removal more efficient and effective.

Our revolutionary Visionaire® Touchscreen Control Interface is integrated into the cabinet to control & monitor all aspects of product operation.



**CARBON
FILTER**

Size Specification (W x D x H)

CT 800	CT 1100	CT 1400	CT 1800
800mm x 700mm x 1285mm	1100mm x 700mm x 1285mm	1400mm x 700mm x 1285mm	1800mm x 700mm x 1285mm
CT PRO 800	CT PRO 1100	CT PRO 1400	CT PRO 1800
800mm x 700mm x 1345mm *	1100mm x 700mm x 1345mm *	1400mm x 700mm x 1345mm *	1800mm x 700mm x 1345mm *

*When Sliding Sash is at Maximum Opening.